

## Instruments Set

**Cat. Ref.: SET204 - Description: Emergency gynecology set**

Items included	Qt.
Sterilizing Stainless Steel Box 30x12x6 cm	1
Metzemaum Scissors Straight 18 cm	1
Operating Scissors B/B Curved 16 cm	1
Kocher Forceps 1x2 Teeth Straight 16 cm	3
Mair Polypus Forceps Straight 24 cm	1
J.L Faure Forceps Curved 1x2 Teeth 22 cm	1
Museux Forceps Straight 24 cm 7 x 5 mm	2
Pozzi Forceps Straight 24 cm	2
Collin Vaginal Speculum 35 mm x 110 mm	1
Cusco Vaginal Speculum pliable 30 mm x 90 mm	1
Probe With Tread 30 cm Uterin	1
Doyen Vaginal Speculum 90 mm	1